NWM038-2 PATIENT TROLLEY





NWM038-2

FEATURES AND HIGHLIGHTS

Ajuste de placa trasera

Technical parameters

External size (LxW)	1960x700 mm
 Table height adjustment 	620-860 mm
Back plate adjustment	0-75°

Technical configuration

🗹 Crank system	1 set
🗹 Back plate air spring	1 set
🗹 Plastic side guardrail	1 set
☑ 5" covered castors	4 pcs
🗹 IV pole	1 set
☑ IV pole prevision	4 pcs
Mattress	1 pc



Guardrail

Virgin ABS material tuck away side rail with buffer. The guardrail is shaped on easy structure to hold and it can help the patient stand up easily by holding the rail.




Handle crank

Steel hanging crank, 100,000 times remaining intact, it can be fold under the bed when not in use.



Bed surface The trolley surface is made of powder coated steel, smooth surface no welding scar.




Bed frame The trolley base is made of powder coated steel.



Caster

5 inch TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts.

Pass Dynamic test: Bearing 120kg run 30KM and pass obstacles 500time.



Steels | Painting Back rest lifting 0-75°, controlled by gas spring.



NUMBER WIN



+86 18021239272

♥ No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

export@number-win.com

() www.number-win.com





As the circumstance of printing, the photos might have minor difference on color, tactile from real products. As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.

General catalogue